

Product datasheet for **RC200969A1V**

Human CPLX1 (NM_006651) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human CPLX1 (NM_006651) AAV Particle

Tag: Myc-DDK

Symbol: CPLX1

Synonyms: CPX-I; CPX1; DEE63; EIEE63

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC200969 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGTTTGTGATGAAGCAGGCTCTAGGAGGGGCCACCAAGGACATGGGAAGATGCTGGGGGGTGACG
AGGAGAAGGACCCAGACGCCGCCAAGAAGGAGGAGGAGCGGCAGGAGGCGCTGCGCCAGGCGGAGGAGGA
GCGCAAGGCCAAGTACGCCAAGATGGAGGCGGAGCGCGAGGCCGTGCGCCAGGGCATCCGAGACAAGTAC
GGCATCAAGAAGAAGGAGGAGCGCGAGGCCGAGGCCAGGCCGCCATGGAGGCCAACTCCGAGGGGAGCT
TGACGCGGCCCAAGAAGGCCATCCCGCCGGGCTGCGGGGACGAGGTGGAGGAGGAGGACGAGAGCATCCT
GGACACCGTCATCAAGTACCTGCCCGGGCCGCTGCAGGACATGCTCAAGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200969 protein sequence
Red=Cloning site Green=Tags(s)

MEFVMKQALGGATKDMGKMLGGDEEKDPDAAKKEEERQEALRQAEERKAKYAKMEAEREAVRQGIRDKY
GIKKKEEREAEAQAAMEANSEGLTRPKKAIPP GCGDEVEEEDSILDTVIKYLPGPLQDMLKK

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_006651



[View online »](#)

ORF Size:	402 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006651.1</u>
RefSeq Size:	2200 bp
RefSeq ORF:	405 bp
Locus ID:	10815
UniProt ID:	<u>O14810</u>
Cytogenetics:	4p16.3
MW:	15 kDa