

Product datasheet for **RC200975L4V**

PDXK (NM_003681) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PDXK (NM_003681) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PDXK
Synonyms:	C21orf97; C21orf124; HEL-S-1a; HMSN6C; PKH; PNK; PRED79
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003681
ORF Size:	936 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200975).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003681.3
RefSeq Size:	7390 bp
RefSeq ORF:	939 bp
Locus ID:	8566
UniProt ID:	O00764
Cytogenetics:	21q22.3
Domains:	pfkB
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Metabolic pathways, Vitamin B6 metabolism

MW: 35.1 kDa

Gene Summary: The protein encoded by this gene phosphorylates vitamin B6, a step required for the conversion of vitamin B6 to pyridoxal-5-phosphate, an important cofactor in intermediary metabolism. The encoded protein is cytoplasmic and probably acts as a homodimer. Alternatively spliced transcript variants have been described, but their biological validity has not been determined. [provided by RefSeq, Jul 2008]