

Product datasheet for **RC200979A1V**

Human JAB1 (COPS5) (NM_006837) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human JAB1 (COPS5) (NM_006837) AAV Particle
Tag:	Myc-DDK
Symbol:	JAB1
Synonyms:	CSN5; JAB1; MOV-34; SGN5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200979 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGGCGTCCGGGAGCGGTATGGCCAGAAAACCTGGGAACGGCCAACAACATGCAGGAAGCTCAGA
GTATCGATGAAATCTACAAATACGACAAGAAACAGCAGCAAGAAATCCTGGCGGCGAAGCCCTGGACTAA
GGATCACCATTACTTTAAGTACTGCAAAATCTCAGCATTGGCTCTGCTGAAGATGGTATGCATGCCAGA
TCGGGAGGCAACTTGAAGTGATGGGTCTGATGCTAGGAAAGGTGGATGGTGAACCATGATCATTATGG
ACAGTTTTGCTTTGCCTGTGGAGGGCACTGAAACCCGAGTAAATGCTCAGGCTGCTGCATATGAATACAT
GGCTGCATACATAGAAAATGCAAAACAGGTTGGCCGCCTTGAAAATGCAATCGGGTGGTATCATAGCCAC
CCTGGCTATGGCTGCTGGCTTTCTGGGATTGATGTTAGTACTCAGATGCTCAATCAGCAGTTCCAGGAAC
CATTTGTAGCAGTGGTGATTGATCCAACAAGAACAATATCCGCAGGGAAGTGAATCTTGGCGCCTTTAG
GACATACCCAAAGGGCTACAAACCTCCTGATGAAGACCTTCTGAGTACCAGACTATTCCACTTAATAAA
ATAGAAGATTTTGGTGACTGCAACAATATTATGCCTTAGAAGTCTCATATTTCAAATCCTCTTTGG
ATCGCAAATTGCTTGAGCTGTTGTGAATAAACTACTGGGTGAATACGTTGAGTTCTTCTAGCTTGCTTAC
TAATGCAGACTATACCACTGGTCAGGTCTTTGATTTGTCTGAAAAGTTAGAGCAGTCAGAAGCCAGCTG
GGACGAGGGAGTTTCATGTTGGGTTTAGAAACGCATGACCGAAAATCAGAAGACAACTTGCCAAAGCTA
CAAGAGACAGCTGTAAACTACCATAGAAGCTATCCATGGATTGATGTCTCAGGTTATTAAGGATAAACT
GTTTAATCAAATTAACATCTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200979 protein sequence Red=Cloning site Green=Tags(s) MAASGSGMAQKTWELANNMQEAQSIDEIYKYDKKQQQEILAAKPWTKDHHYFKYCKISALALLKMMVHAR SGGNLEVMGLMLGKVDGETMIIMDSFALPVEGTETRVNAQAAAYEYMAAYIENAKQVGRLENAIGWYHSH PGYGCWLSGIDVSTQMLNQQFQEPFVAVVIDPRTTISAGKVNLAGFRTYPKGYKPPDEGPSEYQTIPLNK IEDFGVHCKQYYALEVSYFKSSLDKRLLELLWNKYWVNTLSSSSLLTNADYTTGQVFDLSEKLEQSEAQL GRGSFMLGLETHDRKSEDKLAKATRDSCKTTIEAIHGLMSQVIKDKLFNQINIS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006837
ORF Size:	1002 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006837.2</u>
RefSeq Size:	1510 bp
RefSeq ORF:	1005 bp
Locus ID:	10987
UniProt ID:	<u>Q92905</u>
Cytogenetics:	8q13.1
MW:	37.6 kDa