

Product datasheet for RC200980L4

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD71 (TFRC) (NM_003234) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: CD71 (TFRC) (NM_003234) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: CD71

Synonyms: CD71; IMD46; p90; T9; TFR; TFR1; TR; TRFR

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC200980).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_003234

ORF Size: 2280 bp





CD71 (TFRC) (NM_003234) Human Tagged Lenti ORF Clone - RC200980L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 003234.1</u>

 RefSeq Size:
 5010 bp

 RefSeq ORF:
 2283 bp

 Locus ID:
 7037

 UniProt ID:
 P02786

Cytogenetics: 3q29

Domains: PA, TFR_dimer, Peptidase_M28

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protease, Secreted Protein, Transmembrane

Protein Pathways: Endocytosis, Hematopoietic cell lineage

MW: 84.7 kDa


Gene Summary: This gene encodes a cell surface receptor necessary for cellular iron uptake by the process of

receptor-mediated endocytosis. This receptor is required for erythropoiesis and neurologic development. Multiple alternatively spliced variants have been identified. [provided by

RefSeq, Sep 2015]



Product images:



Circular map for RC200980L4