

Product datasheet for **RC200993A1V**

Human MRI1 (NM_001031727) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human MRI1 (NM_001031727) AAV Particle
Tag:	Myc-DDK
Symbol:	MRI1
Synonyms:	M1Pi; MRDI; MTNA; Ypr118w
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200993 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACCCTGGAGGCGATCCGCTACTCGCGGGGCTCCCTGCAGATCCTAGACCAGCTGCTGCTGCCAAGC
AGAGCCGCTACGAGGCGGTGGGCTCGGTGCACCAAGGCTGGGAGGCCATCCGCGCCATGAAGGTGCGGGG
CGCCCCGGCCATAGCCCTGGTGGGCTGTCTCAGCCTCGCCGTGGAGCTGCAGGCGGGCGCCGGGGGACCG
GGACTCGCCGCGCTCGTGGCCTTCGTGCGGACAAGCTGAGCTTCCTCGTCACCGCCCGGCCACCGCTG
TCAACATGGCCCGCGCCCGCCGACCTGGCTGATGTTGCAGCCCGGAGGCCGAACGGGAGGGCGCTAC
GGAAGAGGCGGTCCGGGAGAGAGTGATCTGCTGCACCGAGGACATGCTGGAGAAAGACCTCAGAGACAAC
CGAAGCATTGGGACCTAGGAGCCCGCCACCTCCTGGAGCGGGTGGCCCCAGCGGTGGCAAGGTGACTG
TGCTGACCCACTGTAACACTGGTGTCTGGCCACCGCTGGCTATGGTACAGCCCTAGGTGTGATTCGCTC
ACTGCACAGCCTGGGCCGCTGGAGCATGCCTTCTGCACAGAGACCCGGCCCTACAACCAGGGAGCCCGG
CTGACGGCCTTTGAGCTGGTCTATGAGCAGATCCCCGCCACCCTATCACCGACAGCATGGTGGCTGCTG
CCATGGCCCATAGGGCGGTGTCAGCTGTGGTCTGGGAGCTGACCGCGTGGTTGCCAACGGCGACACAGC
CAACAAGGTGGGCACCTACCAGCTGGCCATTGTCGCCAAGCACCATTGGCATTCCCTTCTACGTGGCTGCC
CCCAGCTCTTCATGTGACCTCCGTCTGGAGACCGCAAGGAGATCATTATTGAAGAGCGACCGGGCCAGG
AGCTGACCGATGTTAATGGGGTCCGGATTGCAGCACCTGGGATTGGAGTTTGGAATCCTGCCTTCGATGT
CACCCCCACGACCTCATCACTGGTGGCATCATCAGAACTGGGGGTCTTTGCCCTGAGGAGCTCCGG
ACAGCCCTAACCAACCATCTCTCCAGGGATGGAACCTAGATGGACCCAGATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



View online »

Protein Sequence:	>RC200993 protein sequence Red=Cloning site Green=Tags(s) MTLEAIRYSRGSLLQILDQLLLKPQSRYEAVGVSQHQAWEAIRAMKVRGAPAIALVGCLSLAVELQAGAGGP GLAALVAFVRDKLSFLVTARPTAVNMAARAARDLADVAAREAEREGATEEAVRERVICCTEDMLEKDLRDN RSIGDLGARHLLERVAPSGGKVTVLTHCNTGALATAGYGTALGVIRSLHSLGRLEHAFCCTETRPYNQGAR LTAFELVYEQIPATLITDSMVAAMAHRGVSAVVVGADRVVANGDTANKVGTYQLAIVAKHHGIPFYVAA PSSSCDLRLETGKEIIIEERPGQELTDVNGVRIAAPGIGVWNPAPFDVTPHDLITGGIITELGVFAPEELR TALTTTISSRDGTLDGPM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001031727
ORF Size:	1107 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001031727.2
RefSeq Size:	3177 bp
RefSeq ORF:	1110 bp
Locus ID:	84245
UniProt ID:	Q9BV20
Cytogenetics:	19p13.13
MW:	39.1 kDa