

## Product datasheet for **RC201004L3V**

### DCTN5 (NM\_032486) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DCTN5 (NM_032486) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DCTN5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032486
ORF Size:	546 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201004).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032486.2</a>
RefSeq Size:	7329 bp
RefSeq ORF:	549 bp
Locus ID:	84516
UniProt ID:	<a href="#">Q9BTE1</a>
Cytogenetics:	16p12.2
Domains:	hexapep
MW:	20.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a subunit of dynactin, a component of the cytoplasmic dynein motor machinery involved in minus-end-directed transport. The encoded protein is a component of the pointed-end subcomplex and is thought to bind membranous cargo. A pseudogene of this gene is located on the long arm of chromosome 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2011]