

## Product datasheet for **RC201007L3V**

### COLEC11 (NM\_024027) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	COLEC11 (NM_024027) Human Tagged ORF Clone Lentiviral Particle
Symbol:	COLEC11
Synonyms:	3MC2; CL-11; CL-K1-I; CL-K1-II; CL-K1-IIa; CL-K1-IIb; CLK1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024027
ORF Size:	813 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201007).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_024027.3</a>
RefSeq Size:	1705 bp
RefSeq ORF:	816 bp
Locus ID:	78989
UniProt ID:	<a href="#">Q9BWP8</a>
Cytogenetics:	2p25.3
MW:	28.7 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the collectin family of C-type lectins that possess collagen-like sequences and carbohydrate recognition domains. Collectins are secreted proteins that play important roles in the innate immune system by binding to carbohydrate antigens on microorganisms, facilitating their recognition and removal. The encoded protein binds to multiple sugars with a preference for fucose and mannose. Mutations in this gene are a cause of 3MC syndrome-2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]