

Product datasheet for **RC201031L3V**

TAOK3 (NM_016281) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TAOK3 (NM_016281) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TAOK3
Synonyms:	DPK; hKFC-A; JIK; MAP3K18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_016281
ORF Size:	2694 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201031).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016281.2
RefSeq Size:	4399 bp
RefSeq ORF:	2697 bp
Locus ID:	51347
UniProt ID:	Q9H2K8
Cytogenetics:	12q24.23
Domains:	pkinese, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

Protein Pathways: MAPK signaling pathway

MW: 105.4 kDa

Gene Summary: The protein encoded by this gene is a serine/threonine protein kinase that activates the p38/MAPK14 stress-activated MAPK cascade but inhibits the basal activity of the MAPK8/JNK cascade. The encoded protein is a member of the GCK subfamily of STE20-like kinases. [provided by RefSeq, Oct 2016]