

Product datasheet for **RC201033A1V**

Human **C11ORF67 (AAMDC) (NM_024684) AAV Particle**

Product data:

Product Type: AAV Particles
Tag: Myc-DDK
Symbol: C11ORF67
Synonyms: C11orf67; CK067; PTD015
Mammalian Cell None
Selection:
Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC201033 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGACTTCCCCTGAAATTGCTTCCTTATCATGGGGCAAATGAAAGTAAAAGGCTCTAATACAACCTATA
AGGACTGCAAAGTATGGCCAGGGGTAGTCGACTTGGGATTGGAGAGAAACAGGAAGTGAACCTCTCC
TGGTGTGCAGCCTGCAGATGTGAAGGAAGTTGTTGAGAAGGGTGTACAGACTCTTGTGATTGGCCGAGGG
ATGAGTGAGGCCTTGAAGGTGCCTTCATCAACTGTGGAGTACCTCAAGAAACATGGCATTGATGTGCGGG
TCCTCCAGACAGAGCAGCAGTGAAGGAGTATAATGCCTTGGTTGCCAAGGGGTGAGGGTGGGAGGTGT
CTTCATTCCACCTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201033 protein sequence
Red=Cloning site Green=Tags(s)

MTSPEIASLSWQMVKVGSNTTYKDCKVWPGGSRWDWRETGTEHSPGVQPADVKEVVEKGVQTLVIGRG
MSEALKVPSSSTVEYLKKHGIDVRVLQTEQAVKEYNALVAQGVVRVGVFHSCTC

TRTRPLE**QKLISEEDLAANDILDYKDDDDK**V

Species: Human
Serotype: AAV-2
ACCN: NM_024684



ORF Size:	366 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_024684.2
RefSeq Size:	557 bp
RefSeq ORF:	369 bp
Locus ID:	28971
UniProt ID:	Q9H7C9
Cytogenetics:	11q14.1
MW:	13.3 kDa