

Product datasheet for **RC201033L4V**

C11ORF67 (AAMDC) (NM_024684) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C11ORF67 (AAMDC) (NM_024684) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C11ORF67
Synonyms:	C11orf67; CK067; PTD015
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024684
ORF Size:	366 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201033).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024684.2
RefSeq Size:	557 bp
RefSeq ORF:	369 bp
Locus ID:	28971
UniProt ID:	Q9H7C9
Cytogenetics:	11q14.1
MW:	13.3 kDa



[View online »](#)

Gene Summary:

May play a role in preadipocyte differentiation and adipogenesis.[UniProtKB/Swiss-Prot Function]