

Product datasheet for **RC201061L4V**

CMTM6 (NM_017801) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CMTM6 (NM_017801) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CMTM6
Synonyms:	CKLFSF6; PRO2219
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017801
ORF Size:	549 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201061).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017801.2
RefSeq Size:	3384 bp
RefSeq ORF:	552 bp
Locus ID:	54918
UniProt ID:	Q9NX76
Cytogenetics:	3p22.3
Protein Families:	Transmembrane
MW:	20.4 kDa



[View online »](#)

Gene Summary:

This gene belongs to the chemokine-like factor gene superfamily, a novel family that is similar to the chemokine and transmembrane 4 superfamilies. This gene is one of several chemokine-like factor genes located in a cluster on chromosome 3. This gene is widely expressed in many tissues, but the exact function of the encoded protein is unknown. [provided by RefSeq, Jul 2008]