

Product datasheet for **RC201062A1V**

Human MED18 (NM_017638) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human MED18 (NM_017638) AAV Particle
Tag: Myc-DDK
Symbol: MED18
Synonyms: p28b; SRB5
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201062 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGGCACCTCCAGTCACCATGATGCCTGTCACTGGGGCACCATTAAACATGATGGAGTACCTGTTGC
AGGGAAGTGTTTTAGATCACAGTTTGAAAGCCTCATCCACCGCCTTCGTGGTTGTGTGACAACATGGA
ACCTGAGACTTTCCTTGACCATGAGATGGTATTCCTCCTTAAGGGCCAGCAAGCCAGCCCATTTGTTCTC
AGGGCCCGACGCTCTATGGACAGGGCAGGGGCACCCTGGCATCTGCGCTACCTGGGACAGCCAGAAATGG
GAGACAAGAACCAGCATGCCCTGGTGCGAACTGCGTGGACATTGCCACATCTGAGAACCTCACCGACTT
CTTGATGGAATGGGCTTCGCATGGACCATGAGTTTGTGCTAAGGGACATTTGTTCCGTAAGGGCATC
ATGAAGATTATGGTGTACAAGATTTCCGCATCCTGGTGCCAGGGAACACAGACAGCACTGAGGCCTTGT
CACTCTCTATCTCGTGAATTAAGTGTGGTAGCACCCGCTGGGCAGGACATGGTCTCTGATGACATGAA
GAACTTCGAGAACAGCTAAACCTCTGGTTCACCTAGAGAAAATAGACCCCAAGAGGCTCATG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201062 protein sequence
Red=Cloning site Green=Tags(s)

MEAPPVTMPVTGGTINMMEYLLQGSVLDHSLESLIHRLRGLCDNMPEPETFLDHEMVFLKQGQASPFVL
RARRSMDRAGAPWHLRYLGQPEMGDKNRHALVRNCVDIATSENLTDFLMEMGFRMDHEFVAKGHLFRKGI
MKIMVYKIFRILVPGNTDSTEALSLSYLVELSVVAPAGQDMVSDDMKNF AEQLKPLVHLEKIDPKRLM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_017638
ORF Size:	624 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_017638.1</u>
RefSeq Size:	1869 bp
RefSeq ORF:	627 bp
Locus ID:	54797
UniProt ID:	<u>Q9BUE0</u>
Cytogenetics:	1p35.3
MW:	23.7 kDa