

Product datasheet for **RC201069A1V**

Human CDC42EP3 (NM_006449) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CDC42EP3 (NM_006449) AAV Particle
Tag:	Myc-DDK
Symbol:	CDC42EP3
Synonyms:	BORG2; CEP3; UB1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201069 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCAGCCAAGACCCCAATTTACCTGAAAGCAGCCAATAACAAGAAAGGAAAGAAATTTAACTGAGGG
ACATTCTGTCTCCTGATATGATCAGTCCCCCGCTTGAGACTTTCGCCACACCATCCACATTGGCAAAGA
GGGCCAGCAGATGTCTTTGGAGATATTTCTTTCTTCAAGGGAACACGAGCTTTTACCTGGAAACCAG
GAGAAAGCACACCTGGGCCAGTTCCTGGGCATAATGAGTTCTTCCGGGCCAACAGCACCTCGGACTCTG
TGTTACAGAAACGCCCTCCCCGGTGCTCAAAAATGCCATCTCCCTCCCGACCATTTGGAGGATCCCAAGC
TCTCATGTTGCCCTTATTGTCAACAGTGACATTTAATTCCAAACAGGAGTCCTTCGGGCCAGCAAAGCTG
CCCAGGCTTAGCTGCGAGCCCGTCATGGAGGAAAAAGCTCAGGAGAAAAGCAGTCTGTTGGAGAATGGGA
CAGTCCACCAGGGAGACACCTCGTGGGGCTCCAGCGTTCTGCATCTCAGTCCAGCCAAGGCAGAGACAG
CCACTCCTCCAGCCTGTCCGAACAGTACCCCGACTGGCCAGCCGAGGACATGTTTGACCATCCCACCCCA
TGCGAGCTCATCAAGGAAAGACTAAGTCAGAGGAGTCCCTCTCTGACCTTACAGGTTCCCTCCTCTCCC
TGCAGCTTGATCTTGGGCCCTCACTTTTGGATGAGGTGCTGAATGTAATGGATAAAAAATAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201069 protein sequence Red=Cloning site Green=Tags(s) MPAKTPIYLKAANNKKGKKFLRDILSPDMISPLGDFRHTIHIGKEGQHDVFGDISFLQGNYELLPGNQ EKAHLGQFPGHNEFFRANSTSDSVFTETPSPVLKNAISLPTIGGSQALMLPLLSPVTFNSKQESFGPAKL PRLSCEPVMEEKAQEKSSLLENGTVHQDTSWGSSGSASQSSQGRDSSHSSSLSEQYPDWPAEDMFDHPTP CELIKGTKSEESLSDLTGSLLSLQLDLGPSLLDEVLNVM DKNK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006449
ORF Size:	762 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006449.3</u>
RefSeq Size:	5715 bp
RefSeq ORF:	765 bp
Locus ID:	10602
UniProt ID:	<u>Q9UKI2</u>
Cytogenetics:	2p22.2
MW:	27.7 kDa