

## Product datasheet for **RC201074L4V**

### TRIM44 (NM\_017583) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRIM44 (NM_017583) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRIM44
Synonyms:	AN3; DIPB; HSA249128; MC7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017583
ORF Size:	1032 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201074).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017583.3</a>
RefSeq Size:	4337 bp
RefSeq ORF:	1035 bp
Locus ID:	54765
UniProt ID:	<a href="#">Q96DX7</a>
Cytogenetics:	11p13
Domains:	zf-B_box
MW:	38.5 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, namely a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. [provided by RefSeq, Jul 2008]