

OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for RC201078L3

Acid Phosphatase (ACP1) (NM_007099) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Acid Phosphatase (ACP1) (NM_007099) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	Acid Phosphatase
Synonyms:	HAAP; LMW-PTP; LMWPTP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201078).
Restriction Sites:	Sgfl-Mlul
Cloning Scheme:	
	Cloning sites used for ORF Shuttling:
	Sgf1 ORF Mlu I GCG ATC GC ATG// NNN ACG CGT



ACCN: ORF Size: NM_007099 474 bp



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


	cid Phosphatase (ACP1) (NM_007099) Human Tagged Lenti ORF Clone – RC201078L3
OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <u>custsupport@origene.com</u> or by calling 301.340.3188 option 3 for pricing and delivery.
	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Met	 centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 007099.2</u>
RefSeq Size:	1568 bp
RefSeq ORF:	477 bp
Locus ID:	52
UniProt ID:	<u>P24666</u>
Cytogenetics:	2p25.3
Domains:	LMWPc
Protein Families:	Druggable Genome, Phosphatase, Transmembrane
Protein Pathways:	Adherens junction, Riboflavin metabolism
MW:	18 kDa

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

Scid Phosphatase (ACP1) (NM_007099) Human Tagged Lenti ORF Clone – RC201078L3

Gene Summary:The product of this gene belongs to the phosphotyrosine protein phosphatase family of
proteins. It functions as an acid phosphatase and a protein tyrosine phosphatase by
hydrolyzing protein tyrosine phosphate to protein tyrosine and orthophosphate. This enzyme
also hydrolyzes orthophosphoric monoesters to alcohol and orthophosphate. This gene is
genetically polymorphic, and three common alleles segregating at the corresponding locus
give rise to six phenotypes. Each allele appears to encode at least two electrophoretically
different isozymes, Bf and Bs, which are produced in allele-specific ratios. Multiple
alternatively spliced transcript variants encoding distinct isoforms have been identified for
this gene. [provided by RefSeq, Aug 2008]

Product images:



Circular map for RC201078L3

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US