

Product datasheet for RC201079A1V

Human SSR3 (NM_007107) AAV Particle

Product data:

Product Type: AAV Particles

Tag: Myc-DDK

Symbol: SSR3

Synonyms: TRAPG

Mammalian Cell None

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC201079 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCTCCTAAAGGCAGCTCCAACAGCAGTCTGAGGAGGACCTGCTCCTGCAGGATTTAGCCGCAATC
TCTCGGCCAAGTCTCCGCGCTCTTCTCGAAACGCGTTCATCGTGTCTGCCATCCCCTCTGGTTATA
CTGGCGAATATGGCATATGGATCTTATTCAGTCTGCTGTTTTGTATAGTGTGATGACCCTAGTAAGCACA
TATTTGGTAGCCTTTGCATACAAGAATGTGAAATTTGTTCTCAAGCACAAAGTAGCACAGAAGAGGGAGG
ATGCTGTTTCCAAGAAGTGACTCGAAAATTTCTGAAGCTGATAATAGAAAGATGTCTCGGAAGGAGAA
AGATGAAAGAATCTTGTGGAAGAAGAATGAAGTTGCTGATTATGAAGCTACAACATTTCCATCTTCTAT
AACAACACTCTGTTCTCGTGGTCTGTCATTGTTGCTTCTTTCATATTGAAGAAGTCAACCCACAG
TGAACACATATTGTCCATAAGTGCTTCATCAGGACTCATCGCCCTCCTGTCTACTGGCTCCAAA

AGCGGACCGACGCGTACGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201079 protein sequence
Red=Cloning site Green=Tags(s)

MAPKGSSKQSEEDLLLQDFSRNLSAKSSALFFGNAFIVSAIPIWLYWRIWHMDLIQSAVLVSVMTLVST
YLVAFAYNKVKFVLKHKVAQKREDAVSKEVTRKLSEADNRKMSRKEKDERILWKKNEVADYEATTSIFY
NNTLFLVVIVASFFILKNFNPTVNYILSISASSGLIALSTGSK

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2



ACCN:	NM_007107
ORF Size:	555 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_007107.2</u>
RefSeq Size:	3734 bp
RefSeq ORF:	558 bp
Locus ID:	6747
UniProt ID:	<u>Q9UNL2</u>
Cytogenetics:	3q25.31
MW:	21.1 kDa