

Human CCDC85B (NM_006848) AAV Particle

Product Type:	AAV Particles
----------------------	---------------

Tag: Myc-DDK

Symbol: CCDC85B

Synonyms: DIPA

Mammalian Cell	None
----------------	------

Selection:

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC201086 representing NM_006848
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGAGGCCGAGGCAGGCGGCCTGGAGGAGCTGACGGACGAGGAGATGGCGGCGCTAGGCAAGGAAGAGC
TAGTGGCGCGCTGCGGCGGGAGGAGGCGGCGCGCTGGCGGCACTGGTGCAGCGCGGCCGCTCATGCA
GGAGTGTAATCGGCAGCTGACAGGGCCACCTGGGCGAGATCCGCGAGCTCAAGCAGCTCAACCGCGCTCTG
CAGGCAGAGAACCCTGAGCTGCGCGACCTCTGCTGCTTCTGGACTCGGAGCGCCAGCGCGGGCGGCGCG
CCGCACGCCAGTGGCAGCTCTTCGGGACCCAAGCATCCCGGGCGTGC GCGAGGACCTGGGCGGCTGTTG
GCAGAAGCTGGCCGAGCTGGAGGGCGGCCAGGAGGAGCTGCTGCGGGAGAACCTAGCGCTTAAGGAGCTC
TGCTGGCGCTGGGCGAAGAATGGGCGTCCCGCGGCGGCCACCGCGCGCCGGGGATCAGGAGCCGGGC
CAGCACCCAGAGCTGCCCTTGCCCGCTGCGGCGCCCGCCAGCTAGGCGATGGAAGCTCCAGCACTGGCAG
CGTGGGCACTCCGAGTCAGTTGCCCTGGCCTGTTCCTCCCGATGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201086 representing NM_006848
Red=Cloning site Green=Tags(s)

MEAEAGGLEELTDEEMAALGKEELVRRLRREEAARLAALVQRGRLMQEVNRQLQGHLEIRELKQLNRRL
QAENRELRDLCCFLDSEQRGRRAARQWQLFGTQASRAVREDLGGCWQKLAELEGRQEELLRENALKELE
CLALGEEWGRPGGSGAGGSGAGAPAPALPPCGPRDLGDGSSSTGSGVSPDOLPLACSPDD

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Species: Human



Serotype:	AAV-2
ACCN:	NM_006848
ORF Size:	606 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006848.2
RefSeq Size:	1277 bp
RefSeq ORF:	609 bp
Locus ID:	11007
UniProt ID:	Q15834
Cytogenetics:	11q13.1
MW:	21.9 kDa