

Product datasheet for **RC201091A1V**

Human **PARM1 (PARM1) (NM_015393) AAV Particle**

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	PARM1
Synonyms:	Cipar1; DKFZP564O0823; PARM-1; WSC4
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC201091 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTCTACAAGACTCTCTCGCTCTTTGCATCTTAACCGCAGGATGGAGGTACAGAGTCTGCCTACAT
CAGCTCCTTTGTCTGTTTCTCTCCGACAAACATTGTACCACCGACCACCATCTGGACTAGCTCTCCACA
AAACACTGATGCAGACTGCCTCCCATCCAACGGCACTCACAACAACCTCGGTGCTCCAGTTACAGCA
TCAGCCCCAACATCTCTGCTTCTTAAGAACATTTCATAGAGTCCAGAGAAGAGGAGATCACCAGCCAG
GTTGCAATTGGGAAGGCACAAACACAGACCCCTCACCTTCTGGGTTCTCGTCAACAAGCGGTGGAGTCCA
CTTAACAACCACGTTGAAGGAACACAGCTCGGGCACTCCTGAAGCAGGCGTGGCAGCTACACTGTCGAG
TCCGCTGCTGAGCCTCCACACTCATCTCCCTCAAGCTCCAGCCTCATCACCTCATCCCTATCAACCT
CACCACCTGAGGTCTTTCTGCCTCCGTTACTACCAACCATAGCTCCACTGTGACCAGCACCCAACCCAC
TGGAGCTCCAACGCACCAGAGTCCCGACAGAGGAGTCCAGCTCTGACCACACCCCACTTCACATGCC
ACAGCTGAGCCAGTACCCAGGAGAAAACACCCCAACAACGTGTGTCAGGCAAGTGATGTGTGAGCTCA
TAGACATGGAGACCACCACCTTTCCAGGGTGATCATGCAGGAAGTAGAATGCATTAAAGTTCAGG
CAGCATCGCCGCCATTACCGTGACAGTCATTGCCGTGGTGTGCTGGTGTGAGTTGCAGCCTACCTA
AAAATCAGGCATTCTCTATGGAAGACTTTTGGACGACCATGACTACGGGTCTGGGAAACTACAACA
ACCCTCTGTACGATGACTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:	>RC201091 protein sequence Red=Cloning site Green=Tags(s) MYYKTLFALCILTAGWRVQSLPTSAPLSVSLPTNIVPPTTIWTSSPQNTDADTASPSNGTHNNSVLPVTA SAPTSLLPKNISIESREEEITSPGSNWEGTNDPSPSGFSSTSGGVHLTTTLKEHSSGTPEAGVAATLSQ SAAEPPTLISPQAPASSPSSLSTSPPEVFSASVTTNHSSTVTSTQPTGAPTAPESPTEESSDHTPTSHA TAEPVPQEKTPPTTVSGKVMCELDIMETTTTFPRVIMQEVEHALSSGSIAAITVTVIAVLLVFGVAAYL KIRHSSYGRLLDDHDYGSWGNYNPLYDDS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015393
ORF Size:	930 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_015393.2
RefSeq Size:	5069 bp
RefSeq ORF:	933 bp
Locus ID:	25849
UniProt ID:	Q6UWI2
Cytogenetics:	4q13.3
MW:	32.3 kDa