

Product datasheet for RC201094

PECI (ECI2) (NM_206836) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PECI (ECI2) (NM_206836) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PECI
Synonyms:	ACBD2; dj1013A10.3; DRS-1; DRS1; HCA88; PECI
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC201094 representing NM_206836 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGATGGCGTACTTGGCTTGGAGACTGGCGGGCGTTCTGTCCGAGTTCTCTGCAGGTCAGT
TCCCGGTAGTTCAGCTGCACATGAATAGAACAGCAATGAGAGCCAGTCAGAAGGACTTTGAAAATCAAT
GAATCAAGTGAACTCTTAAAAAGGATCCAGGAAACGAAGTGAAGCTAAAACCTACGCGCTATATAAG
CAGGCCACTGAAGGACCTTGTAACTGCCCAAACAGGTGATTTGACTTGATCAACAAGGCCAAATGGG
ACGCATGGAATGCCCTTGGCAGCCTGCCAAGGAAGCTGCCAGGCAGAACTATGTGGATTTGGTGTCCAG
TTTGAGTCCTTCAATTGGAATCCTCTAGTCAGGTGGAGCCTGGAACAGACAGGAAATCAACTGGGTTTGAA
ACTCTGGTGGTGACCTCCGAAGATGGCATCACAAAGATCATGTTCAACCGGCCCAAAAAGAAAAATGCCA
TAAACTGAGATGTATCATGAAATTAATGCGTGCACCTAAAGCTGCCAGCAAGGATGACTCAATCATCAC
TGTTTTAACAGGAAATGGTGACTATTACAGTAGTGGGAATGATCTGACTAACTTCACTGATATCCCCCT
GGTGGAGTAGAGGAGAAAGCTAAAAATAATGCCGTTTTACTGAGGGAATTTGTGGGCTGTTTTATAGATT
TTCCTAAGCCTCTGATTGCAGTGGTCAATGGTCCAGCTGTGGGCATCTCCGTCACCCTCCTGGGCTATT
CGATGCCGTGTATGCATCTGACAGGGCAACATTTACATACACCATTTAGTCACCTAGGCCAAAGTCCGGAA
GGATGCTCCTCTTACACTTTTCCGAAGATAATGAGCCCAGCCAAGGCAACAGAGATGCTTATTTTTGGAA
AGAAGTTAACAGCGGGAGAGGCATGTGCTCAAGGACTTGTACTGAAGTTTTCCCTGATAGCACTTTTCA
GAAAGAAGTCTGACCAGGCTGAAGGCATTTGCAAAGCTTCCCCCAAATGCCTTGAGAATTTCAAAGAG
GTAATCAGGAAAAGAGAGAGAGAAAAACTACACGCTGTTAATGCTGAAGAATGCAATGTCCTTCAGGGAA
GATGGCTATCAGATGAATGCACAAATGCTGTGGTGAACCTCTTATCCAGAAAATCAAACCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC201094 representing NM_206836
 Red=Cloning site Green=Tags(s)

MAMAYLAWRLARRSCPSSLQVTSFPVVQLHMNRTAMRASQKDFENSMNQVKLLKKDPGNEVKLKLALYK
 QATEGPCNMPKPGVFDL INKAKWDANWALGSLPKAARQNYVDL VSSLSPLSESSQVPEGTDRKSTGFE
 TLVVTSEDGITKIMFNRPKKKNAINTEMYHEIMRALKAASKDSSIITVL TGNGDYSSGNDLTNFTDIPP
 GGVEEKAKNNAVLLREFVGCIDFPKPLIAVNGPAVGISVTLGLFDVAVYASDRATFHTPF SHLGGQSPE
 GCSSYTFPKIMSPAKATEMLIFGKLTAGEACAQGLVTEVFPDSTFQKEVWTRLKAFKLPNALRISKE
 VIRKREREKLVAVNAEECNVLQGRWLSDECTNAVNVNFLSRKSKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_206836

ORF Size: 1182 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_206836.3](#)

RefSeq Size: 1420 bp

RefSeq ORF: 1185 bp

Locus ID: 10455

UniProt ID: [O75521](#)


Cytogenetics: 6p25.2

Protein Pathways: Fatty acid metabolism

MW: 43.6 kDa

Gene Summary: This gene encodes a member of the hydratase/isomerase superfamily. The protein encoded is a key mitochondrial enzyme involved in beta-oxidation of unsaturated fatty acids. It catalyzes the transformation of 3-cis and 3-trans-enoyl-CoA esters arising during the stepwise degradation of cis-, mono-, and polyunsaturated fatty acids to the 2-trans-enoyl-CoA intermediates. Alternatively spliced transcript variants have been described. [provided by RefSeq, Aug 2011]

Product images:



Circular map for RC201094