

Product datasheet for RC201095

MDH2 (NM_005918) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MDH2 (NM_005918) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MDH2
Synonyms:	DEE51; EIEE51; M-MDH; MDH; MGC:3559; MOR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC201095 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTCTCCGCCCTCGCCGGCCTGTCAGCGCTGCTCTCCGCCGAGCTTCAGCACCTCGGCCAGAACA
ATGCTAAAGTAGCTGTGCTAGGGCCTCTGGAGGCATCGGGCAGCCACTTCACTTCTCTGAAGAACAG
CCCCTTGGTGAGCCGCTGACCTCTATGATATCGCGCACACCCGGAGTGGCCGAGATCTGAGCCAC
ATCGAGACCAAAGCCGCTGTAAAGGCTACCTCGGACCTGAACAGCTGCCTGACTGCCTGAAAGGTTGTG
ATGTGGTAGTTATTCGGCTGGAGTCCCAGAAAGCCAGGCATGACCCGGGACGACCTGTCAACACCAA
TGCCACGATTGTGGCCACCCTGACCGCTGCCTGTGCCAGCACTGCCCGAAGCCATGATCTGCCTCATT
GCCAATCCGGTTAATTCCACCATCCCCATCACAGCAGAAGTTTTCAAGAAGCATGGAGTGTACAACCCCA
ACAAAACTTTCGGCGTGACGACCCTGGACATCGTCAGAGCCAACACCTTTGTTGCAGAGCTGAAGGGTTT
GGATCCAGCTCGAGTCAACGTCCCTGTCATTGGTGGCCATGCTGGGAAGACCATCATCCCCCTGATCTCT
CAGTGCACCCCAAGGTGGACTTCCCCAGGACCAGCTGACAGCACTCACTGGGCGGATCCAGGAGGCCG
GCACGGAGGTGGTCAAGGCTAAAGCCGGAGCAGGCTCTGCCACCCTCTCCATGGCGTATGCCGGCAGCCG
CTTTGTCTTCTCCCTTGTGGATGCAATGAATGAAAGGAAGGTGTTGTGGAATGTTCTTCGTTAAGTCA
CAGGAAACGGAATGTACCTACTTCTCACACCGCTGCTGCTTGGGAAAAAGGGCATCGAGAAGAACCTGG
GCATCGGCAAAGTCTCCTTTTGAGGAGAAGATGATCTCGGATGCCATCCCCGAGCTGAAGGCCTCCAT
CAAGAAGGGGAAGATTTCTGTAAGACCCTGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC201095 protein sequence
Red=Cloning site Green=Tags(s)

MLSALARPVSAALRRSFSTSAQNNKAVLVGASGGIGQPLSLLLKNSPLVSRLTLYDIAHTPGVAADLSH
 IETKAAVKGYLGPEQLPDCLKGCDVVVIPAGVPRKPGMTRDDL FNTNATIVATLTAACAQHCPEAMICVI
 ANPVNSTIPITAEVFKKHGVYNPKNKIFGVTTLDIVRANTFVAELKGLDPARVNPVIGGHAGKTIIP LIS
 QCTPKVDFPQDQLTALTGRIQEAGTEVVKAKAGAGSATLSMAYAGARFVSLVDAMNGKEGVVECSFVKS
 QTECTYFSTPLLLGKKGIEKNLIGIKVSSFEEKMISDAIPELKASIKKGEDFVKTLK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk6083_b10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




Plasmid Map:




ACCN: NM_005918

ORF Size:	1014 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_005918.4
RefSeq Size:	2268 bp
RefSeq ORF:	1017 bp
Locus ID:	4191
UniProt ID:	P40926
Cytogenetics:	7q11.23
Domains:	ldh
Protein Families:	Druggable Genome
Protein Pathways:	Citrate cycle (TCA cycle), Glyoxylate and dicarboxylate metabolism, Metabolic pathways, Pyruvate metabolism
MW:	35.5 kDa
Gene Summary:	Malate dehydrogenase catalyzes the reversible oxidation of malate to oxaloacetate, utilizing the NAD/NADH cofactor system in the citric acid cycle. The protein encoded by this gene is localized to the mitochondria and may play pivotal roles in the malate-aspartate shuttle that operates in the metabolic coordination between cytosol and mitochondria. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2013]

Product images:



Western blot validation of overexpression lysate (Cat# [LY401789]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC201095 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified MDH2 protein (Cat# [TP301095]). The protein was produced from HEK293T cells transfected with MDH2 cDNA clone (Cat# RC201095) using MegaTran 2.0 (Cat# [TT210002]).