

## Product datasheet for **RC201098L1V**

### SF3A1 (NM\_005877) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SF3A1 (NM_005877) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SF3A1
Synonyms:	PRP21; PRPF21; SAP114; SF3A120
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005877
ORF Size:	2379 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201098).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005877.4</a>
RefSeq Size:	5159 bp
RefSeq ORF:	2382 bp
Locus ID:	10291
UniProt ID:	<a href="#">Q15459</a>
Cytogenetics:	22q12.2
Domains:	Surp, UBP
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Spliceosome

**MW:** 88.9 kDa

**Gene Summary:** This gene encodes a subunit of the splicing factor 3a protein complex. The splicing factor 3a heterotrimer is a component of the mature U2 small nuclear ribonucleoprotein particle (snRNP). U2 small nuclear ribonucleoproteins play a critical role in spliceosome assembly and pre-mRNA splicing. [provided by RefSeq, Aug 2014]