

Product datasheet for **RC201117A1V**

Human TFIIA2 (GTF2A2) (NM_004492) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TFIIA2 (GTF2A2) (NM_004492) AAV Particle
Tag: Myc-DDK
Symbol: TFIIA2
Synonyms: HsT18745; T18745; TF2A2; TFIIA; TFIIA-12; TFIIA-gamma; TFIAS
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201117 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCATATCAGTTATACAGAAATACTACTTTGGGAAACAGTCTTCAGGAGAGCCTAGATGAGCTCATAC
 AGTCTCAACAGATCACCCCCCACTTGCCCTTCAAGTTCTACTTCAGTTTGATAAGGCTATAATGCAGC
 ACTGGCTCAGAGGGTCAGGAACAGAGTCAATTCAGGGGCTCTCTAAATACGTACAGATTCTGCGATAAT
 GTGTGGACTTTTGTACTGAATGATGTTGAATTCAGAGAGGTGACAGAATTATTAAAGTGGATAAAGTGA
 AAATTGTAGCCTGTGATGGTAAAAATACTGGCTCCAATACTACAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201117 protein sequence
 Red=Cloning site Green=Tags(s)
 MAYQLYRNTTLGNSLQESLDELISQSQITPQLALQVLLQFDKAINAALAQRVRNRVNFGRSLNTYRFCDN
 VWTFLVNDVEFREVTELIKVDKVKIVACDGKNTGSNTTE

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_004492
ORF Size: 327 bp


[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004492.1
RefSeq Size:	1582 bp
RefSeq ORF:	330 bp
Locus ID:	2958
UniProt ID:	P52657
Cytogenetics:	15q22.2
MW:	12.5 kDa