

Product datasheet for RC201119A1V

Human DUSP3 (NM_004090) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human DUSP3 (NM_004090) AAV Particle
Tag: Myc-DDK
Symbol: DUSP3
Synonyms: VHR
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201119 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGTCGGGCTCGTTTCGAGCTCTCGGTGCAGGATCTCAACGACCTGCTCTCGGACGGCAGCGGCTGCTACA
 GCCTCCCGAGCCAGCCCTGCAACGAGGTACCCCGCGGATCTACGTGGGCAACGCGTCTGTGGCTCAGGA
 CATCCCAAGCTGCAGAACTAGGCATACCCATGTGCTGAACGCGGCTGAGGGCAGGTCCTTCATGCAC
 GTCAACACCAATGCCAACTTCTACAAGGACTCCGGCATCACATACCTGGGCATCAAGGCCAACGACACAC
 AGGAGTTCAACCTCAGCGTTACTTTGAAAGGGCTGCCGACTTCATTGACCAGGCTTTGGCTCAAAAGAA
 TGGCCGGGTGCTCGTCCACTGCCGGAAGGTTATAGCCGCTCCCAACGCTAGTTATCGCTACCTCATG
 ATGCGGCAGAAGATGGACGTCAAGTCTGCCCTGAGCATCGTGAGGCAGAACCGTGAGATCGGCCCAACG
 ATGGCTTCTGGCCAGCTCTGCCAGCTCAATGACAGACTAGCCAAGGAGGGGAAGTTGAAACCC

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201119 protein sequence
 Red=Cloning site Green=Tags(s)

MSGSFELSVQDLNDLLSDGSGCYSLPSQPCNEVTPRIYVGNASVAQDIPKLQKLGIHVLNAAEGRSFMH
 VNTNANFYKDSGITYLGIKANDTQEFNLSAYFERAADFIDQALAQKNGRVLVHCREGYSRSPTLVIAYLM
 MRQKMDVKSALSIVRQNREIGPNDGFLAQLCQLNDRLAKEGKLKP

SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_004090
ORF Size:	555 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004090.2
RefSeq Size:	4139 bp
RefSeq ORF:	558 bp
Locus ID:	1845
UniProt ID:	P51452
Cytogenetics:	17q21.31
MW:	20.5 kDa