

Product datasheet for **RC201120L3V**

DAP1 (DAP) (NM_004394) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DAP1 (DAP) (NM_004394) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DAP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004394
ORF Size:	306 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201120).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004394.1
RefSeq Size:	2385 bp
RefSeq ORF:	309 bp
Locus ID:	1611
UniProt ID:	P51397
Cytogenetics:	5p15.2
MW:	11.2 kDa


[View online »](#)

Gene Summary:

This gene encodes a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death that is induced by interferon-gamma. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, May 2014]