

## Product datasheet for **RC201127L2V**

### **KEL (NM\_000420) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KEL (NM_000420) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KEL
Synonyms:	CD238; ECE3; Kell
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_000420
ORF Size:	2196 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201127).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000420.2</a>
RefSeq Size:	2562 bp
RefSeq ORF:	2199 bp
Locus ID:	3792
UniProt ID:	<a href="#">P23276</a>
Cytogenetics:	7q34
Domains:	Peptidase_M13
Protein Families:	Druggable Genome, Protease, Transmembrane



[View online »](#)

**MW:** 82.8 kDa

**Gene Summary:** This gene encodes a type II transmembrane glycoprotein that is the highly polymorphic Kell blood group antigen. The Kell glycoprotein links via a single disulfide bond to the XK membrane protein that carries the Kx antigen. The encoded protein contains sequence and structural similarity to members of the neprilysin (M13) family of zinc endopeptidases. [provided by RefSeq, Jul 2008]