

## **Product datasheet for RC201132L2**

#### 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com

OriGene Technologies, Inc.

EU: info-de@origene.com CN: techsupport@origene.cn

## UQCRC2 (NM\_003366) Human Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

Product Name: UQCRC2 (NM 003366) Human Tagged Lenti ORF Clone

Tag: mGFP
Symbol: UQCRC2

Synonyms: MC3DN5; QCR2; UQCR2

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

**E. coli Selection:** Chloramphenicol (34 ug/mL)

**ORF Nucleotide** 


The ORF insert of this clone is exactly the same as(RC201132).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





st The last codon before the Stop codon of the ORF.

**ACCN:** NM\_003366

ORF Size: 1359 bp





### UQCRC2 (NM\_003366) Human Tagged Lenti ORF Clone - RC201132L2

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

> reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: NM 003366.2

RefSeq Size: 1674 bp RefSeq ORF: 1362 bp Locus ID: 7385 **UniProt ID:** P22695

Cytogenetics:

16p12.2 Domains: Peptidase M16, Peptidase M16 C

**Protein Families:** Druggable Genome, Protease

**Protein Pathways:** Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways,

Oxidative phosphorylation, Parkinson's disease

MW: 48.4 kDa


**Gene Summary:** The protein encoded by this gene is located in the mitochondrion, where it is part of the

> ubiquinol-cytochrome c reductase complex (also known as complex III). This complex constitutes a part of the mitochondrial respiratory chain. Defects in this gene are a cause of

mitochondrial complex III deficiency nuclear type 5. [provided by RefSeq, Jul 2015]



# **Product images:**



Circular map for RC201132L2