

## Product datasheet for **RC201132L3V**

### UQCRC2 (NM\_003366) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	UQCRC2 (NM_003366) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UQCRC2
Synonyms:	MC3DN5; QCR2; UQCR2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003366
ORF Size:	1359 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201132).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003366.2</a>
RefSeq Size:	1674 bp
RefSeq ORF:	1362 bp
Locus ID:	7385
UniProt ID:	<a href="#">P22695</a>
Cytogenetics:	16p12.2
Domains:	Peptidase_M16, Peptidase_M16_C
Protein Families:	Druggable Genome, Protease



[View online »](#)

**Protein Pathways:** Alzheimer's disease, Cardiac muscle contraction, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

**MW:** 48.4 kDa

**Gene Summary:** The protein encoded by this gene is located in the mitochondrion, where it is part of the ubiquinol-cytochrome c reductase complex (also known as complex III). This complex constitutes a part of the mitochondrial respiratory chain. Defects in this gene are a cause of mitochondrial complex III deficiency nuclear type 5. [provided by RefSeq, Jul 2015]