

Product datasheet for **RC201136A1V**

Human CRALBP (RLBP1) (NM_000326) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CRALBP (RLBP1) (NM_000326) AAV Particle
Tag:	Myc-DDK
Symbol:	CRALBP
Synonyms:	CRALBP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201136 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCAGAAGGGGTGGGCACGTTCCGCATGGTACCTGAAGAGGAACAGGAGCTCCGTGCCCAACTGGAGC
AGCTCACAAACCAAGGACCATGGACCTGTCTTTGGCCCGTGACGCCAGCTGCCCCGCCACACCTTGCAGAA
GGCCAAGGATGAGCTGAACGAGAGAGAGGAGACCCGGGAGGAGGCAGTGCGAGAGCTGCAGGAGATGGTG
CAGGCGCAGGCGGCCTCGGGGGAGGAGCTGGCGGTGGCCGTGGCGGAGAGGGTGCAAGAGAAGGACAGCG
GCTTCTTCTCGCTTCATCCGCGCACGGAAGTTCAACGTGGGCCGTGCCTATGAGCTGCTCAGAGGCTA
TGTGAATTTCCGGCTGCAGTACCCTGAGCTCTTTGACAGCCTGTCCCCAGAGGCTGTCCGCTGCACCAT
GAAGCTGGCTACCCTGGTGTCTCTCTAGTCGGGACAAGTATGGCCGAGTGGTCATGCTCTTCAACATTG
AGAAGCTGGCAAAGTCAAGAAATCACCTTTGATGAGATCTTGACGGCATATTGCTTCATCCTGGAGAAGCT
GCTGGAGAATGAGGAACTCAAATCAATGGCTTCTGCATCATTGAGAATTCAAGGGCTTTACCATGCAG
CAGGCTGCTAGTCTCCGACTTCAGATCTCAGGAAGATGGTGGACATGCTCCAGGATTCTTCCCAGCCC
GGTTCAAAGCCATCCACTTCATCCACCAGCCATGGTACTTCACCACGACCTACAATGTGGTCAAGCCCTT
CTTGAAGAGCAAGCTGCTTGAGAGGGTCTTTGTCCACGGGGATGACCTTTCTGGTTTCTACCAGGAGATC
GATGAGAACATCCTGCCCTCTGACTTCGGGGGCACGCTGCCCAAGTATGATGGCAAGGCCGTTGCTGAGC
AGCTCTTTGGCCCCCAGGCCCAAGCTGAGAACACAGCCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201136 protein sequence Red=Cloning site Green=Tags(s) MSEGVGTFRMPVEEQELRAQLEQLTTKDHGPVFGPCSQLPRHTLQKAKDELNEREETREEAVRELQEMV QAQAASGEELAVAVAERVQEKDSGFFLRFI RARKFNVGRAYELLRGYVNFRLQYPELFDSL SPEAVRCTI EAGYPGVLSSRDKYGRVVMLFNIENWQSQEITFDEILQAYCFILEKLLENEETQINGFCIIENFKGFTMQ QAASLRTSDLRKMVDMLQDSFPARFKAIHFIIHQPWYFTTTYNVVKPFLKSKLLERVFVHGDDL SGFYQEI DENILPSDFGGTLPKYDGKAVAEQLFGPQAQAENTAF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_000326
ORF Size:	951 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_000326.3</u>
RefSeq Size:	1752 bp
RefSeq ORF:	954 bp
Locus ID:	6017
UniProt ID:	<u>P12271</u>
Cytogenetics:	15q26.1
MW:	36.5 kDa