

Product datasheet for **RC201138A1V**

Human RFC2 (NM_002914) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RFC2 (NM_002914) AAV Particle
Tag:	Myc-DDK
Symbol:	RFC2
Synonyms:	RFC40
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201138 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGGTGGAGGCCGTCTGTGGTGGCGCGGGCGAGGTGGAGGCCAGGACTCTGACCCTGCCCTGCCT
TCAGCAAGGCCCCCGGCAGCGCCGGCCACTACGAACTGCCGTGGGTGAAAAATATAGGCCAGTAAAGCT
GAATGAAATTGTCGGGAATGAAGACACCGTGAGCAGGCTAGAGGTCTTTGCAAGGAAGGAAATGTGCC
AACATCATCATTGCGGGCCCTCCAGGAACCGGCAAGACCACAAGCATTCTGTGCTTGGCCCGGCCCTGC
TGGGCCCAGCACTCAAAGATGCCATGTTGGAACCTCAATGCTTCAAATGACAGCATGACCGACGGAGCCCA
GCAAGCCTTGAGGAGAACCATGGAAATCTACTCTAAAACCACTCGCTTCGCCCTTGCTTGAATGCTTCG
GATAAGATCATCGAGCCATTCACTCCCGCTGTGAGTCCTCCGGTACACAAAGCTGACCGACGCCGAGA
TCCTCACCAGGCTGATGAATGTTATCGAGAAGGAGAGGGTACCCTACCTGATGACGGCCTAGAAGCCAT
CATCTTCACGGCCAGGGAGACATGAGGCAGGCGCTGAACAACCTGCAGTCCACCTTCTCAGGATTTGGC
TTCATTAACAGTGAGAACGTGTTCAAGGTCTGTGACGAGCCCCACCACTGCTGGTAAAGGAGATGATCC
AGCACTGTGTGAATGCCAACATTGACGAAGCCTACAAGATTCTTGCTCACTTGTGGCATCTGGGCTACTC
ACCAGAAGATATCATTGGCAACATCTTCGAGTGTGTAACCTTTCAAATGGCAGAATACCTGAAACTG
GAGTTTATCAAGGAAATTGGATACACTCACATGAAATAGCGGAAGGAGTGAACCTCTTTTGCAGATGG
CAGGCCTCTGGCAAGGCTGTGTGAGAACAATGGCCCCGGTGGCCAGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201138 protein sequence Red=Cloning site Green=Tags(s) MEVEAVCGGAGEVEAQSDPAPAFSKAPGSAGHYELPWVEKYRPVKLNEIVGNEDTVSRLEVFAREGNVP NIIIAGPPGTGKTTTILCLARALLGPALKDAMLELNASNDSDGAAQALRRRTMEIYSKTTTFALACNAS DKIIEPIQSRCAVLRYSKLTDAQILTRLMNVIEKERVPTDDGLEAIIIFTAQGDMRQALNNLQSTFSGFG FINSENVFKVCDEPHLLVKEMIQHCVNANIDEAYKILAHLWHLGYSPEDIIGNIFRVCKTFQMAEYLLK EFIKEIGYTHMKIAEGVNSLLQAGLLARLCQKTMAPVAS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002914
ORF Size:	960 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_002914.3
RefSeq Size:	1657 bp
RefSeq ORF:	963 bp
Locus ID:	5982
UniProt ID:	P35250
Cytogenetics:	7q11.23
MW:	35.2 kDa