

Product datasheet for **RC201138L4V**

RFC2 (NM_002914) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RFC2 (NM_002914) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RFC2
Synonyms:	RFC40
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002914
ORF Size:	960 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201138).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002914.3
RefSeq Size:	1657 bp
RefSeq ORF:	963 bp
Locus ID:	5982
UniProt ID:	P35250
Cytogenetics:	7q11.23
Domains:	AAA, AAA
Protein Families:	Druggable Genome, Stem cell - Pluripotency



[View online »](#)

Protein Pathways: DNA replication, Mismatch repair, Nucleotide excision repair

MW: 35.2 kDa

Gene Summary: This gene encodes a member of the activator 1 small subunits family. The elongation of primed DNA templates by DNA polymerase delta and epsilon requires the action of the accessory proteins, proliferating cell nuclear antigen (PCNA) and replication factor C (RFC). Replication factor C, also called activator 1, is a protein complex consisting of five distinct subunits. This gene encodes the 40 kD subunit, which has been shown to be responsible for binding ATP and may help promote cell survival. Disruption of this gene is associated with Williams syndrome. Alternatively spliced transcript variants encoding distinct isoforms have been described. A pseudogene of this gene has been defined on chromosome 2. [provided by RefSeq, Jul 2013]