

Product datasheet for RC201152

IDH2 (NM_002168) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: IDH2 (NM_002168) Human Tagged ORF Clone

Tag: Myc-DDK

Symbol: IDH2

Synonyms: D2HGA2; ICD-M; IDH; IDHM; IDP; IDPM; mNADP-IDH

Mammalian Cell Neomycin

Selection:

Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn





ORF Nucleotide Sequence:

>RC201152 ORF sequence

Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

CGGCCCTGACAGCCCCCACCTCGCAAGAGCAGCCGCGGCGCCACTATGCCGACAAAAGGATCAAGGTGGC GAAGCCCGTGGTGGAGATGGATGATGAGATGACCCGTATTATCTGGCAGTTCATCAAGGAGAAGCTC ATCCTGCCCCACGTGGACATCCAGCTAAAGTATTTTGACCTCGGGCTCCCAAACCGTGACCAGACTGATG ACCAGGTCACCATTGACTCTGCACTGGCCACCCAGAAGTACAGTGTGGCTGTCAAGTGTGCCACCATCAC CCCTGATGAGGCCCGTGTGGAAGAGTTCAAGCTGAAGAAGATGTGGAAAAGTCCCAATGGAACTATCCGG AACATCCTGGGGGGGACTGTCTTCCGGGAGCCCATCATCTGCAAAAACATCCCACGCCTAGTCCCTGGCT GGACCAAGCCCATCACCATTGGCAGGCACGCCCATGGCGACCAGTACAAGGCCACAGACTTTGTGGCAGA CCGGGCCGGCACTTTCAAAATGGTCTTCACCCCAAAAGATGGCAGTGGTGTCAAGGAGTGGGAAGTGTAC AACTTCCCCGCAGGCGCGTGGGCATGGGCATGTACAACACCGACGAGTCCATCTCAGGTTTTGCGCACA GCTGCTTCCAGTATGCCATCCAGAAGAAATGGCCGCTGTACATGAGCACCAAGAACACCATACTGAAAGC CTACGATGGCCGTTTCAAGGACATCTTCCAGGAGATCTTTGACAAGCACTATAAGACCGACTTCGACAAG AATAAGATCTGGTATGAGCACCGGCTCATTGATGACATGGTGGCTCAGGTCCTCAAGTCTTCGGGTGGCT TTGTGTGGGCCTGCAAGAACTATGACGGAGATGTGCAGTCAGACATCCTGGCCCAGGGCTTTGGCTCCCT TGGCCTGATGACGTCCGTCCTGGTCTGCCCTGATGGGAAGACGATTGAGGCTGAGGCCGCTCATGGGACC GTCACCCGCCACTATCGGGAGCACCAGAAGGGCCGGCCCACCAGCACCAACCCCATCGCCAGCATCTTTG CCTGGACACGTGGCCTGGAGCACCGGGGGAAGCTGGATGGGAACCAAGACCTCATCAGGTTTGCCCAGAT GCTGGAGAAGGTGTGCGTGGAGACGGTGGAGAGTGGAGCCATGACCAAGGACCTGGCGGGCTGCATTCAC GGCCTCAGCAATGTGAAGCTGAACGAGCACTTCCTGAACACCACGGACTTCCTCGACACCATCAAGAGCA ACCTGGACAGAGCCCTGGGCAGGCAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC201152 protein sequence
Red=Cloning site Green=Tags(s)

MAGYLRVVRSLCRASGSRPAWAPAALTAPTSQEQPRRHYADKRIKVAKPVVEMDGDEMTRIIWQFIKEKL ILPHVDIQLKYFDLGLPNRDQTDDQVTIDSALATQKYSVAVKCATITPDEARVEEFKLKKMWKSPNGTIR NILGGTVFREPIICKNIPRLVPGWTKPITIGRHAHGDQYKATDFVADRAGTFKMVFTPKDGSGVKEWEVY NFPAGGVGMGMYNTDESISGFAHSCFQYAIQKKWPLYMSTKNTILKAYDGRFKDIFQEIFDKHYKTDFDK NKIWYEHRLIDDMVAQVLKSSGGFVWACKNYDGDVQSDILAQGFGSLGLMTSVLVCPDGKTIEAEAAHGT VTRHYREHQKGRPTSTNPIASIFAWTRGLEHRGKLDGNQDLIRFAQMLEKVCVETVESGAMTKDLAGCIH GLSNVKLNEHFLNTTDFLDTIKSNLDRALGRQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk6153 a02.zip

Restriction Sites: Sgfl-Mlul



Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_002168

ORF Size: 1356 bp

OTI Disclaimer:

Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at customport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <u>More info</u>

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
- 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.



RefSeq: <u>NM 002168.4</u>

 RefSeq Size:
 1818 bp

 RefSeq ORF:
 1359 bp

 Locus ID:
 3418

 UniProt ID:
 P48735

 Cytogenetics:
 15q26.1

 Domains:
 isodh

Protein Pathways: Citrate cycle (TCA cycle), Glutathione metabolism, Metabolic pathways

MW: 50.9 kDa


Gene Summary: Isocitrate dehydrogenases catalyze the oxidative decarboxylation of isocitrate to 2-

oxoglutarate. These enzymes belong to two distinct subclasses, one of which utilizes NAD(+) as the electron acceptor and the other NADP(+). Five isocitrate dehydrogenases have been reported: three NAD(+)-dependent isocitrate dehydrogenases, which localize to the mitochondrial matrix, and two NADP(+)-dependent isocitrate dehydrogenases, one of which is mitochondrial and the other predominantly cytosolic. Each NADP(+)-dependent isozyme is a homodimer. The protein encoded by this gene is the NADP(+)-dependent isocitrate dehydrogenase found in the mitochondria. It plays a role in intermediary metabolism and

energy production. This protein may tightly associate or interact with the pyruvate dehydrogenase complex. Alternative splicing results in multiple transcript variants. [provided


by RefSeq, Feb 2014]


Product images:



Circular map for RC201152







Western blot validation of overexpression lysate (Cat# [LY400787]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC201152 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).

Coomassie blue staining of purified IDH2 protein (Cat# [TP301152]). The protein was produced from HEK293T cells transfected with IDH2 cDNA clone (Cat# RC201152) using MegaTran 2.0 (Cat# [TT210002]).