

Product datasheet for **RC201182L4V**

LDB1 (NM_003893) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LDB1 (NM_003893) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LDB1
Synonyms:	CLIM-2; CLIM2; LDB-1; NLI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_003893
ORF Size:	1125 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201182).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003893.3
RefSeq Size:	2305 bp
RefSeq ORF:	1128 bp
Locus ID:	8861
UniProt ID:	Q86U70
Cytogenetics:	10q24.32
Domains:	LIM_bind
Protein Families:	Transcription Factors



[View online »](#)

MW: 42.8 kDa

Gene Summary: Binds to the LIM domain of a wide variety of LIM domain-containing transcription factors. May regulate the transcriptional activity of LIM-containing proteins by determining specific partner interactions. Plays a role in the development of interneurons and motor neurons in cooperation with LHX3 and ISL1. Acts synergistically with LHX1/LIM1 in axis formation and activation of gene expression. Acts with LMO2 in the regulation of red blood cell development, maintaining erythroid precursors in an immature state (By similarity). [UniProtKB/Swiss-Prot Function]