

## Product datasheet for **RC201185L3V**

### TIPRL (NM\_001031800) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | TIPRL (NM_001031800) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | TIPRL  |
| Synonyms:                 | TIP; TIP41; TIPRL1   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_001031800   |
| ORF Size:                 | 534 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC201185).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_001031800.1</a>   |
| RefSeq Size:              | 954 bp   |
| RefSeq ORF:               | 537 bp   |
| Locus ID:                 | 261726   |
| UniProt ID:               | <a href="#">O75663</a>   |
| Cytogenetics:             | 1q24.2   |
| MW:                       | 20.2 kDa   |



[View online »](#)

**Gene Summary:**

TIPRL is an inhibitory regulator of protein phosphatase-2A (PP2A) (see PPP2CA; MIM 176915), PP4 (see PPP4C; MIM 602035), and PP6 (see PPP6C; MIM 612725) (McConnell et al., 2007 [PubMed 17384681]).[supplied by OMIM, Nov 2010]