

Product datasheet for **RC201201L1V**

PSMD11 (NM_002815) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PSMD11 (NM_002815) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PSMD11
Synonyms:	p44.5; Rpn6; S9
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_002815
ORF Size:	1266 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201201).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002815.2
RefSeq Size:	1598 bp
RefSeq ORF:	1269 bp
Locus ID:	5717
UniProt ID:	O00231
Cytogenetics:	17q11.2
Domains:	PCI
Protein Families:	Stem cell - Pluripotency




[View online »](#)

Protein Pathways: Proteasome

MW: 47.3 kDa

Gene Summary: The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. This gene encodes a member of the proteasome subunit S9 family that functions as a non-ATPase subunit of the 19S regulator and is phosphorylated by AMP-activated protein kinase. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012]

Product images:



[RC201201L1] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC201201L1V particle to overexpress human PSMD11-Myc-DDK fusion protein.