

## Product datasheet for **RC201207L3V**

### CCDC106 (NM\_013301) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CCDC106 (NM_013301) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CCDC106
Synonyms:	HSU79303; ZNF581
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013301
ORF Size:	840 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201207).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013301.2</a> , <a href="#">NP_037433.2</a>
RefSeq Size:	1993 bp
RefSeq ORF:	843 bp
Locus ID:	29903
UniProt ID:	<a href="#">Q9BWC9</a>
Cytogenetics:	19q13.42
MW:	32 kDa



[View online »](#)

**Gene Summary:** Promotes the degradation of p53/TP53 protein and inhibits its transactivity.[UniProtKB/Swiss-Prot Function]