

Product datasheet for **RC201210A1V**

Human TRIB2 (NM_021643) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TRIB2 (NM_021643) AAV Particle
Tag:	Myc-DDK
Symbol:	TRIB2
Synonyms:	C5FW; GS3955; TRB2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201210 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAACATACACAGGTCTACCCCATCACAATAGCGAGATATGGGAGATCGCGGAACAAACCCAGGATT
TCGAAGAGTTGTCGTCTATAAGGTCCGCGGAGCCAGCCAGAGTTTCAGCCGAACCTCGGCTCCCCGAG
CCCGCCGAGACTCCGAACCTGTCGCATTGCGTTTCTTGATCGGGAATACTTATTGTTGGAACCTCTG
GAGGGAGACCACGTTTTTCGTGCCGTGCATCTGCACAGCGGAGAGGAGCTGGTGTGAAGGTGTTTGATA
TCAGCTGCTACCAGGAATCCCTGGCACCCTGCTTTTGCCTGTCTGCTCATAGTAACATCAACCAATCAC
TGAAATTATCCTGGGTGAGACCAAGCCTATGTGTTCTTTGAGCGAAGCTATGGGGACATGCATTCCCTTC
GTCCGCACCTGCAAGAAGCTGAGAGAGGAGGAGGCAGCCAGACTGTTCTACCAGATTGCCTCGGCAGTGG
CCCACTGCCATGACGGGGGGCTGGTGTGCGGGACCTCAAGCTGCGGAAATTCATCTTTAAGGACGAAGA
GAGGACTCGGGTCAAGCTGGAAGCCTGGAAGACGCCTACATTCTCGGGGAGATGATGATTCCCTCTCC
GACAAGCATGGCTGCCCGGCTTACGTAAGCCAGAGATCTTGAACACCAGTGGCAGCTACTCGGGCAAAG
CAGCCGACGTGTGGAGCCTGGGGGTGATGCTGTACACCATGTTGGTGGGCGGTACCTTTCCATGACAT
TGAACCCAGCTCCCTCTTCAGCAAGATCCGGCGTGGCCAGTTCAACATTCAGAGACTCTGTGCCCCAAG
GCCAAGTGCCTCATCCGAAGCATTCTGCGTCGGGAGCCCTCAGAGCGGCTGACCTCGCAGGAAATTTCTGG
ACCATCCTTGGTTTTCTACAGATTTTAGCGTCTCGAATTCAGCATATGGTGCTAAGGAAGTGTCTGACCA
GCTGGTGCCGACGTCAACATGGAAGAGAACTTGACCTTTCTTTAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201210 protein sequence Red=Cloning site Green=Tags(s) MNIHRSTPITIARYGRSRNKTQDFEELSSIRSAEPSQSFSPNLGSPSPPETPNLSHCVSCIGKYLLLEPL EGDHVFRVHLHSGEELVCKVFDISCYQESLAPCFCLSAHSNINQITEIILGETKAYVFFERSYGDMHSF VRTCKKLREEEAARLFYQIASAVAHCHDGGVLRLDLKLRKFIKDEERTRVKLESLEDAYILRGDDDSL DKHGCPAYVSP EILNTSGSYSGKAADVWSLGVMLYTMLVGRYPFHDI EPSSLFSKIRRGQFNIPETLSPK AKCLIRSILRREP SERLTSQEILDHPWFSTDFSVSNSAYGAKEVSDQLVPDVNMEENLDPFFN TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_021643
ORF Size:	1029 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_021643.1
RefSeq Size:	4408 bp
RefSeq ORF:	1032 bp
Locus ID:	28951
UniProt ID:	Q92519
Cytogenetics:	2p24.3
MW:	38.8 kDa