

Product datasheet for **RC201217A1V**

Human CRIP1 (NM_001311) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CRIP1 (NM_001311) AAV Particle
Tag: Myc-DDK
Symbol: CRIP1
Synonyms: CRHP; CRIP; CRP-1; CRP1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201217 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCCCAAGTGTCCCAAGTGAACAAGGAGGTGTACTTCGCCGAGAGGGTGACCTCTCTGGGCAAGGACT
GGCATCGGCCCTGCCTGAAGTGCAGAAATGTGGGAAGACGCTGACCTCTGGGGCCACGCTGAGCACGA
AGGCAAACCCTACTGAACCACCCTGCTACGCAGCCATGTTTGGGCCTAAAGGCTTTGGGCGGGCGGA
GCCGAGAGCCACACTTTCAAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201217 protein sequence
Red=Cloning site Green=Tags(s)

MPKCPKCNKEVYFAERVTSLGKDWHRPCLKCEKCGKLTSGGHAEHEGKPYCNHPCYAAMFGPKGFGRG
AESHTFK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_001311
ORF Size: 231 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001311.3</u>
RefSeq Size:	480 bp
RefSeq ORF:	234 bp
Locus ID:	1396
UniProt ID:	<u>P50238</u>
Cytogenetics:	14q32.33
MW:	8.5 kDa