

Product datasheet for **RC201222A1V**

Human PLAUR (NM_002659) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PLAUR (NM_002659) AAV Particle
Tag:	Myc-DDK
Symbol:	PLAUR
Synonyms:	CD87; U-PAR; UPAR; URKR
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201222 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGTCACCCGCCGCTGCTGCCGCTGCTGCTGCTGCTCCACACCTGCGTCCCAGCCTCTTGGGGCCTGC
GGTGATGTCAGTGTAAAGACCAACGGGGATTGCCGTGTGGAAGAGTGCGCCCTGGGACAGGACCTCTGCAG
GACCACGATCGTGCCTTGTGGGAAGAAGGAGAAGAGCTGGAGCTGGTGGAGAAAAGCTGTACCCACTCA
GAGAAGACCAACAGGACCCTGAGCTATCGGACTGGCTTGAAGATCACCAGCCTTACCGAGGTTGTGTGTG
GGTTAGACTTGTGCAACAGGGCAACTCTGGCCGGGCTGTACCTATTCGGAAGCCGTTACCTCGAATG
CATTTCTGTGGCTCATCAGACATGAGCTGTGAGAGGGGCCGGCACCAGAGCCTGCAGTGCCGCAGCCCT
GAAGAACAGTGCCTGGATGTGGTGACCCACTGGATCCAGGAAGGTGAAGAAGGGCGTCCAAGGATGACC
GCCACCTCCGTGGCTGTGGCTACCTTCCCGGCTGCCCGGGCTCCAATGGTTTCCACAACAACGACACCTT
CCACTTCTGAAATGCTGCAACACCACCAAAATGCAACGAGGGCCCAATCCTGGAGCTTGAAAATCTGCCG
CAGAATGGCCGCCAGTGTACAGCTGCAAGGGGAACAGCACCCATGGATGCTCCTCTGAAGAGACTTTCC
TCATTGACTGCCGAGGCCCCATGAATCAATGTCTGGTAGCCACCGGCACTCACGAACCGAAAAACCAAG
CTATATGGTAAGAGGCTGTGCAACCGCCTCAATGTGCCAATGCCCACCTGGGTGACGCCTTCAGCATG
AACCACATTGATGTCTCTGCTGTACTAAAAGTGGCTGTAACCAACCCAGACCTGGATGTCCAGTACCGCA
GTGGGGCTGCTCCTCAGCCTGGCCCTGCCATCTCAGCCTCACCATCACCTGCTAATGACTGCCAGACT
GTGGGGAGGCACTCTCCTCTGGACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201222 protein sequence Red=Cloning site Green=Tags(s) MGHPPLLPLLLLHTCPASWGLRCMQCKTNGDCRVEECALGQDLCRTTIVRLWEEGEELELVEKSCTHS EKTNRTLSYRTGLKITSLTEVVCGLDLNQGNSGRAVTYSRSRYLECISGSSDMSCERGRHQLQCRSP EEQCLDVVTHWIQEGEEGRPKDDRHLRGCGYLPGCPGSGNFHNDTFHFLKCCNTTKCNEGPILELENLP QNGRQCYSCKGNSTHGCSSEETFLIDCRGPMNQCLVATGTHEPKNQSYMVRGCATASMCQHAHLGDAFSM NHIDVSCCTKSGCNHPDLVDQYRSGAAPQPGPAHLSLTITLLMTARLWGGTLLWT TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002659
ORF Size:	1005 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_002659.2
RefSeq Size:	1570 bp
RefSeq ORF:	1008 bp
Locus ID:	5329
UniProt ID:	Q03405
Cytogenetics:	19q13.31
MW:	37 kDa