

Product datasheet for **RC201227L1V**

SFRS4 (SRSF4) (NM_005626) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SFRS4 (SRSF4) (NM_005626) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SFRS4
Synonyms:	SFRS4; SRP75
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005626
ORF Size:	1482 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201227).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005626.3
RefSeq Size:	2158 bp
RefSeq ORF:	1485 bp
Locus ID:	6429
UniProt ID:	Q08170
Cytogenetics:	1p35.3
Domains:	RRM
Protein Pathways:	Spliceosome



[View online »](#)

MW: 56.5 kDa

Gene Summary: This gene encodes a member of the arginine/serine-rich splicing factor family. The encoded protein likely functions in mRNA processing. [provided by RefSeq, Feb 2009]