

Product datasheet for **RC201230L3V**

PDHB (NM_000925) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PDHB (NM_000925) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PDHB |
| Synonyms: | PDHBD; PDHE1-B; PDHE1B; PHE1B |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_000925 |
| ORF Size: | 1077 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC201230). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_000925.3 , NP_000916.2 |
| RefSeq Size: | 1544 bp |
| RefSeq ORF: | 1080 bp |
| Locus ID: | 5162 |
| UniProt ID: | P11177 |
| Cytogenetics: | 3p14.3 |
| Domains: | transket_pyr, transketolase_C |



[View online »](#)

| | |
|--------------------------|--|
| Protein Pathways: | Butanoate metabolism, Citrate cycle (TCA cycle), Glycolysis / Gluconeogenesis, Metabolic pathways, Pyruvate metabolism, Valine, leucine and isoleucine biosynthesis |
| MW: | 39.2 kDa |
| Gene Summary: | The pyruvate dehydrogenase (PDH) complex is a nuclear-encoded mitochondrial multienzyme complex that catalyzes the overall conversion of pyruvate to acetyl-CoA and carbon dioxide, and provides the primary link between glycolysis and the tricarboxylic acid (TCA) cycle. The PDH complex is composed of multiple copies of three enzymatic components: pyruvate dehydrogenase (E1), dihydrolipoamide acetyltransferase (E2) and lipoamide dehydrogenase (E3). The E1 enzyme is a heterotetramer of two alpha and two beta subunits. This gene encodes the E1 beta subunit. Mutations in this gene are associated with pyruvate dehydrogenase E1-beta deficiency. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Mar 2012] |