

Product datasheet for **RC201231A1V**

Human CD7 (NM_006137) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CD7 (NM_006137) AAV Particle
Tag:	Myc-DDK
Symbol:	CD7
Synonyms:	GP40; LEU-9; Tp40; TP41
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201231 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCGGGCTCCGAGGCTCCTGCTGCTGCCCTGCTTCTGGCGCTGGCTCGCGGCTGCCTGGGGCC
TGGTGGCCAAGAGGTGCAGCAGTCTCCCACTGCACGACTGTCCCGTGGGAGCCTCCGTCAACATCAC
CTGCTCCACGAGCGGGGCTGCGTGGGATCTACCTGAGGCAGCTCGGGCCACAGCCCAAGACATCATT
TACTACGAGGACGGGTGGTGCCCACTACGGACAGACGGTTCGGGGCCGCATCGACTTCTCAGGGTCCC
AGGACAACCTGACTATCACCATGCACCGCTGCAGCTGTGGACACTGGCACCTACACCTGCCAGGCCAT
CACGGAGGTCAATGTCTACGGCTCCGGCACCTGGTCTGTTGACAGAGGAACAGTCCCAAGGATGGCAC
AGATGCTCGGACGCCCCACCAAGGGCTCTGCCCTCCCTGCCCAACGACAGGCTCCGCCCTCCCTGACC
CGCAGACAGCCTCTGCCCTCCCTGACCGCCAGCAGCCTCTGCCCTCCCTGCGGCCCTGGCGGTGATCTC
CTTCTCCTCGGGCTGGGCTGGGGTGGCGTGTGTGCTGGCAGGACACAGATAAAGAACTGTGCTCG
TGGCGGGATAAGAATTCGGCGGCATGTGTGGTGTACGAGGACATGTGCACAGCCGCTGCAACACGCTGT
CCTCCCCCAACCACTACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201231 protein sequence Red=Cloning site Green=Tags(s) MAGPPRLLLLPLLLALARGLPALAAQEVQQSPHCTTVPVGASVNITCSTSGGLRGIYLRQLGPQPQDII YYEDGVVPTTDRRFRGRIDFSGSQDNLITMHRQLSDTGTYTCQAITEVNVYSGTLVLVTEEQSQGWH RCS DAPPRASALPAPPTGSALPDPQTASALPDPPAASALPAALAVISFLLGLGLGVACVLARTQIKKLC WRDKNSAACVVYEDMSHSRCNTLSSPNQYQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006137
ORF Size:	720 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006137.6</u>
RefSeq Size:	1318 bp
RefSeq ORF:	723 bp
Locus ID:	924
UniProt ID:	<u>P09564</u>
Cytogenetics:	17q25.3
MW:	25.4 kDa