

Product datasheet for **RC201231L4V**

CD7 (NM_006137) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | CD7 (NM_006137) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | CD7 |
| Synonyms: | GP40; LEU-9; Tp40; TP41 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_006137 |
| ORF Size: | 720 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC201231). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_006137.6 |
| RefSeq Size: | 1318 bp |
| RefSeq ORF: | 723 bp |
| Locus ID: | 924 |
| UniProt ID: | P09564 |
| Cytogenetics: | 17q25.3 |
| Domains: | ig, IGv, IG |
| Protein Families: | Druggable Genome, Transmembrane |



[View online »](#)

Protein Pathways: Hematopoietic cell lineage

MW: 25.4 kDa

Gene Summary: This gene encodes a transmembrane protein which is a member of the immunoglobulin superfamily. This protein is found on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development. [provided by RefSeq, Jul 2008]