

Product datasheet for RC201238L4

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

POLD1 (NM_002691) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: POLD1 (NM_002691) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: POLD1

Synonyms: CDC2; CRCS10; MDPL; POLD

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC201238).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_002691

ORF Size: 3321 bp





POLD1 (NM_002691) Human Tagged Lenti ORF Clone - RC201238L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 002691.1</u>

 RefSeq Size:
 3464 bp

 RefSeq ORF:
 3324 bp

 Locus ID:
 5424

 UniProt ID:
 P28340

Cytogenetics: 19q13.33

Domains: DNA_pol_B_exo, DNA_pol_B, POLBc

Protein Families: Druggable Genome, Stem cell - Pluripotency

Protein Pathways: Base excision repair, DNA replication, Homologous recombination, Metabolic pathways,

Mismatch repair, Nucleotide excision repair, Purine metabolism, Pyrimidine metabolism

MW: 123.6 kDa


Gene Summary: This gene encodes the 125-kDa catalytic subunit of DNA polymerase delta. DNA polymerase

delta possesses both polymerase and 3' to 5' exonuclease activity and plays a critical role in DNA replication and repair. Alternatively spliced transcript variants have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 6.


[provided by RefSeq, Mar 2012]



Product images:



Circular map for RC201238L4



Double digestion of RC201238L4 using Sgfl and Mlul $\,$