

Product datasheet for **RC201241L2V**

n-Myc (MYCN) (NM_005378) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	n-Myc (MYCN) (NM_005378) Human Tagged ORF Clone Lentiviral Particle
Symbol:	n-Myc
Synonyms:	bHLHe37; MODED; N-myc; NMYC; ODED
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_005378
ORF Size:	1392 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201241).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005378.4
RefSeq Size:	2736 bp
RefSeq ORF:	1395 bp
Locus ID:	4613
UniProt ID:	P04198
Cytogenetics:	2p24.3
Domains:	HLH, Myc_N_term
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 49.6 kDa

Gene Summary: This gene is a member of the MYC family and encodes a protein with a basic helix-loop-helix (bHLH) domain. This protein is located in the nucleus and must dimerize with another bHLH protein in order to bind DNA. Amplification of this gene is associated with a variety of tumors, most notably neuroblastomas. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2014]