

Product datasheet for **RC201250A1V**

Human SNRPD2 (NM_177542) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SNRPD2 (NM_177542) AAV Particle
Tag: Myc-DDK
Symbol: SNRPD2
Synonyms: Sm-D2; SMD2; SNRPD1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201250 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCCTCTCAACAAGCCCAAGAGTGAGATGACCCAGAGGAGCTGCAGAAGCGAGAGGAGGAGGAAT
TTAACACCGGTCCACTCTCTGTGCTCACACAGTCAGTCAAGAACAATACCCAAGTGCTCATCACTGCCG
CAACAATAAGAACTCCTGGGCCGCGTGAAGGCCTTCGATAGGCACTGCAACATGGTGTGGAGAACGTG
AAGGAGATGTGGACTGAGGTACCCAAGAGTGGAAGGGCAAGAAGAAGTCCAAGCCAGTCAACAAGACC
GCTACATCTCCAAGATGTTCTGCGCGGGGACTCAGTCATCGTGGTCTGCGGAACCCGCTCATCGCCGG
CAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201250 protein sequence
 Red=Cloning site Green=Tags(s)

MSLLNPKSEMTPEELQKREEEEFNTGPLSVLTQSVKNNTQVLINCRNNKLLGRVKAFDRHCNMVLENV
KEMWTEVPKSGKGGKSKPVNKDRYISKMFLRGDSVIVVLRNPLIAGK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_177542



[View online »](#)

| | |
|----------------------|---|
| ORF Size: | 354 bp |
| Buffer: | PBS with 0.001% Pluronic F68 |
| Stability: | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| RefSeq: | <u>NM_177542.1</u> , <u>NP_808210.1</u> |
| RefSeq Size: | 882 bp |
| RefSeq ORF: | 327 bp |
| Locus ID: | 6633 |
| UniProt ID: | <u>P62316</u> |
| Cytogenetics: | 19q13.32 |
| MW: | 13.5 kDa |