

## Product datasheet for **RC201273L3V**

### YL1 (VPS72) (NM\_005997) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | YL1 (VPS72) (NM_005997) Human Tagged ORF Clone Lentiviral Particle   |
| Symbol:                   | YL1  |
| Synonyms:                 | CFL1; Swc2; TCFL1; YL-1; YL1   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_005997  |
| ORF Size:                 | 1092 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC201273).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_005997.1</a>  |
| RefSeq Size:              | 1546 bp  |
| RefSeq ORF:               | 1095 bp  |
| Locus ID:                 | 6944   |
| UniProt ID:               | <a href="#">Q15906</a>   |
| Cytogenetics:             | 1q21.3   |
| Protein Families:         | Transcription Factors  |
| MW:                       | 40.6 kDa   |


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a shared subunit of two multi-component complexes, the histone acetyltransferase complex TRRAP/TIP60 as well as the chromatin remodeling SRCAP-containing complex. The TRRAP/TIP60 complex acetylates nucleosomal histones important for transcriptional regulation, double strand DNA break repair and apoptosis. The SRCAP-containing complex catalyzes the exchange of histone H2A with the histone variant Htz1 (H2AFZ) into nucleosomes. This protein may be responsible for binding H2AFZ, which has a role in chromosome segregation. This protein may also have a role in regulating long-term hematopoietic stem cell activity. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Aug 2012]