

## Product datasheet for **RC201283L3V**

### HJURP (NM\_018410) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HJURP (NM_018410) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HJURP
Synonyms:	FAKTS; hFLEG1; URLC9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018410
ORF Size:	2244 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201283).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018410.3</a>
RefSeq Size:	3189 bp
RefSeq ORF:	2247 bp
Locus ID:	55355
UniProt ID:	<a href="#">Q8NCD3</a>
Cytogenetics:	2q37.1
Protein Families:	Druggable Genome
MW:	83.5 kDa



[View online »](#)

**Gene Summary:**

Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at centromeres. Acts as a specific chaperone for CENPA and is required for the incorporation of newly synthesized CENPA molecules into nucleosomes at replicated centromeres. Prevents CENPA-H4 tetramerization and prevents premature DNA binding by the CENPA-H4 tetramer. Directly binds Holliday junctions.[UniProtKB/Swiss-Prot Function]