

Product datasheet for **RC201307L4V**

DHRS11 (NM_024308) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DHRS11 (NM_024308) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DHRS11
Synonyms:	ARPG836; SDR24C1; spDHRS11
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024308
ORF Size:	780 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201307).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024308.3
RefSeq Size:	1608 bp
RefSeq ORF:	783 bp
Locus ID:	79154
UniProt ID:	Q6UWP2
Cytogenetics:	17q12
Domains:	adh_short
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 28.3 kDa

Gene Summary: Catalyzes the conversion of the 17-keto group of estrone, 4- and 5-androstenes and 5-alpha-androstanes into their 17-beta-hydroxyl metabolites and the conversion of the 3-keto group of 3-, 3,17- and 3,20- diketosteroids into their 3-hydroxyl metabolites. Exhibits reductive 3-beta-hydroxysteroid dehydrogenase activity toward 5-beta-androstanes, 5-beta-pregnanes, 4-pregnanes and bile acids. May also reduce endogenous and exogenous alpha-dicarbonyl compounds and xenobiotic alicyclic ketones.[UniProtKB/Swiss-Prot Function]