

Product datasheet for **RC201314A1V**

Human Bcl x (BCL2L1) (NM_138578) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Bcl x (BCL2L1) (NM_138578) AAV Particle
Tag: Myc-DDK
Symbol: Bcl x
Synonyms: Bcl-X; BCL-XL/S; BCL2L; BCLX; PPP1R52
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201314 ORF sequence
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTCTCAGAGCAACCGGGAGCTGGTGGTTGACTTTCTCTCTACAAGCTTTCCAGAAAGGATACAGCT
GGAGTCAGTTTGTGATGTGGAAGAGAACAGGACTGAGGCCCCAGAAGGGACTGAATCGGAGATGGAGAC
CCCCAGTGCCATCAATGGCAACCCATCCTGGCACCTGGCAGACAGCCCCGGTGAATGGAGCCACTGGC
CACAGCAGCAGTTTGGATGCCCGGGAGGTGATCCCCATGGCAGCAGTAAAGCAAGCGCTGAGGGAGGCAG
GCGACGAGTTTGAAGTCCGGTACCGCGGGCATTGACGTGACATCCAGCTCCACATCACCACGAGG
GACAGCATATCAGAGCTTTGAACAGGTAGTGAATGAAGTCTTCCGGGATGGGGTAAAGTGGGGTCCGATT
GTGGCCTTTTTCTCTTCCGGCGGGCACTGTGCGTGGAAAGCGTAGACAAGGAGATGCAGGTATTGGTGA
GTCCGATCGCAGCTTGGATGGCCACTTACCTGAATGACCACCTAGAGCCTTGGATCCAGGAGAACGGCGG
CTGGGATACTTTGTGGAAGTCTATGGGAACAATGCAGCAGCCGAGAGCCGAAAGGGCCAGGAACGCTTC
AACCGCTGGTTCCTGACGGGCATGACTGTGGCCGGCGTGGTTCGCTGGGCTCACTCTTCAGTCGGAAA

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201314 protein sequence
Red=Cloning site **Green**=Tags(s)

MSQSNREL VVDFLSYKLSQKGYSWQFSDVEENRTEAPEGTESEMETPSAINGNPSWHLADSPAVNGATG
HSSSLDAREVIPMAAVKQALREAGDEFELRYRAFSDLTSQLHITPGTAYQSFEQVNVNLFDRDGVNWGRI
VAFFSFGGALCVESVDKEMQVLVSRIAAMATYLNHLEPWIQENGWDTFVELYGNNAAESRKQGERF
NRWFLTGMTVAGVLLGSLFSRK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_138578
ORF Size:	699 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_138578.1
RefSeq Size:	2575 bp
RefSeq ORF:	702 bp
Locus ID:	598
UniProt ID:	Q07817
Cytogenetics:	20q11.21
MW:	26 kDa