

Product datasheet for **RC201316L3V**

C20orf4 (AAR2) (NM_015511) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C20orf4 (AAR2) (NM_015511) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C20orf4
Synonyms:	C20orf4; CGI-23
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015511
ORF Size:	1152 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201316).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015511.2
RefSeq Size:	2765 bp
RefSeq ORF:	1155 bp
Locus ID:	25980
UniProt ID:	Q9Y312
Cytogenetics:	20q11.23
MW:	43.5 kDa



[View online »](#)

Gene Summary:

This gene encodes the homolog of the yeast A1-alpha2 repressin protein that is involved in mRNA splicing. Alternately spliced transcript variants have been found for this gene. [provided by RefSeq, Dec 2012]